Installation Guide for KWedge for Android (7.x)



Version 11/2022

Copyright © 2022 by MSS Software All Rights Reserved

Table of Contents

Overview	2
Steps for Installing and Running KWedge for Android	
Step 1: Install KWedge	
Step 2: Change Folder Permission to Allow Full Control	
Step 3: Select Target	
Step 4: Dock VistA Scanner and Transfer	5
Maximo Instructions	7
Recovering Scanned Data	
Special Messages	
Advanced Options	
Issues with Transferring Data	
Troubleshooting KWedge	
Troubleshooting Checklist	
Scan Count Pro aka SC Pro	
Support	

Overview

KWedge for Android, also referred to as KWedge 7.x, is a data transfer utility for PCs to receive data from barcode scanners provided by MSS Software. Its function allows data to be transferred character by character and line by line, from a mobile barcode reader to a target application on the PC such as Excel or an enterprise application such as VistA. The data will appear in the selected PC application as if it is being typed into the program. KWedge for Android provides an excellent means of transfer where other methods such as WIFI and Bluetooth cannot be used due to security concerns.

KWedge for Android is designed for use with Android based inventory scanners provided by MSS Software and it is not compatible with earlier scanner models such as the CK3X. Models such as the CK3X use a different operating system called Windows Mobile 6.5 and MSS Software provided KWedge 3.23 for use with the Windows based scanners. VA Hospital users are advised to remove KWedge 3.23 from the computer where an Android scanner such as the CT45 Pocket Scanner is used so as to avoid confusion by the user about which KWedge version to use.

KWedge was created by MSS Software and it is designed to be installed on Windows 10 PCs for data communications between the VistA Scanners and the VA's VistA system. Since the VistA Scanners have inventory applications pre-installed, the only thing that needs to be installed by the IT Department is the KWedge for Android software on the Windows 10 PC's. This manual will explain the installation.

The KWedge installation file may be accessed from the MSS Software website or installed from the scanner. Visit the Downloads section of this webpage to find KWedge.



Figure 1: See file path where the KWedge application may be installed from

https://www.msssoftware.com/shop/va_scanning_products/M010AA~CK65+VistA+Scanner+ Kit.htm

The KWedge installation file may also be found on the scanner. The file is located on each scanner from this file path like shown in the screenshot above: <CT45/CK65> \ Internal Shared Storage \ MSS \ KWedgeInstall.exe.

Overview of Steps

- 1. Install KWedge for Android
- 2. Give **Full Permission to delete** data in the install folder (WSS Software\KWedge for Android).
- 3. Start KWedge in Admin mode and select Target (VistA for most users)
- 4. Collect Data on VistA Scanner and select 2. Transfer
- 5. Dock VistA Scanner. On KWedge select Transfer to initiate data transfer

Step 1: Install KWedge

The install file for KWedge for Android is found at the MSS Software website or on the scanner as explained in the previous section. Install on the PC accepting all defaults. Please note the new version of KWedge for Android has a green icon as opposed to the yellow icon (KWedge 3.23) for the older Windows Scanners. Only the Android scanners will work with KWedge for Android. It is recommended that KWedge 3.23 (yellow icon version) be uninstalled from the PC so the user is not confused about which version to use.



Step 2: Change Folder Permission to Allow Full Control

Change the Permission level of the KWedge Installation folder (\MSS Software\KWedge for Android) from Read Only to **Allow Full Control** to delete data. This will allow the transferred file from the scanners (vista.txt.txt or vista.txt) to be deleted after it has been processed.

Step 3: Select Target

The first time KWedge is installed, it is recommended to run the program as administrator so the target application will be pre-selected after each following use. If this step is skipped, the user will be required to Select Target each time KWedge is launched and before data can be transmitted.

Once Run As Administrator has been done as shown in the image below, select the target so that data from the VistA Scanner can be transferred into the VistA or Maximo system by default each time KWedge is launched. Once the reader is docked and a connection has been established, data can easily be transferred to VistA. You do not have to run as administrator from this point on since the target has been established.

Open Move to One Open file loca Run as admin	Drive ation istrator	
(KWedge™ File View		
B	8	
Transfer Dev	ice Time	Action
Select target	10:12a	KWedge for Android version 7.0
	10:12a 10:12a	Current working directory = C:\MSS Software\KWedge for Android. You appear to be operating in a non-administrative mode. To modify options, run as an administrator
3	10:12a	Target application is *Untitled - Notepad
•	10:12a	Target set to 'Excel'

Click the "Select Target" button in KWedge.



The Target window will show any applications that are currently active. If the VistA or Maximo software application is not active prior to launching the KWedge application, then start it VistA or Maximo, and click on the Refresh button in KWedge to make VistA or Maximo visible and select it. (You may also experiment by selecting a program like Excel to test the transfer first).

Example when choosing AbovePar:



Example when choosing MAXIMO:

Select Target	×
Current: MAXIMD New target:	
KWedgel MAXIMO	OK
Microsoft Text Input Application Program Manager	Cancel

Example when choosing Reflections Workspace - VISTA

Select Target	×
Current Reflections Workspace - VISTA New target:	
KWedgel Microsoft Text Input Application Program Manager Reflections Workspace - VISTA	OK Cancel

Step 4: Dock VistA Scanner and Transfer

This step requires that you have an Android VistA Scanner that has collected some data for one of these programs:

- ENNX
- MAXIMO
- PRCUS
- PRCPH
- Scan Count Pro

Once you have collected data on the scanner, go back to the Main Menu and select **2. Transfer Data** which creates an inventory file that can be read by KWedge for Android.

After pressing **2. Transfer Data**, the scanner will ask you to confirm whether or not you want to create a file which can be ready by KWedge. If you choose **Proceed**, it is not possible to review or edit the data on the scanner. Users familiar with the CK3 scanners should take note this is a different process because the CK3 scanners did allow the user to edit their data on the scanner.

Press 2. Transfer Data

Press Proceed when ready to transfer data from scanner to PC

Place scanner in dock. Open KWedge and click Transfer.



Step 5: Transfer Data through KWedge

Make sure the program you are transferring the data into, (VistA, Maximo, Excel or other), is active and ready to receive data. **On KWedge, click the Transfer button**. Data will be transferred to the active program and a pop-up screen showing that the records have been transferred will appear.

🍪 KWedge™	KWedge	\times
File View Transfer Select target	?	VISTA processed (42 records sent; please check results). Clear all transactions from device file now?
		<u>Y</u> es <u>N</u> o

KWedge will then ask if you wish to delete the data from the VistA Scanner. If you are certain that the data has been transferred correctly then click on "Yes" to erase the data. Otherwise, click "No" and then try again.

Maximo Instructions

Steps for Successful Transfer to Maximo

- 1. Launch the SMI (Scanner to Maximo Interface) software program
- 2. Launch the KWedge for Android application
- 3. Within KWedge, choose Select Target, then select SMI (Scanner to Maximo Interface) as shown in Figure 2
- 4. On the scanner, be sure data has already been collection using the Maximo app and Option 2 Transfer Data has been completed so a data file has been created. Place the scanner in its dock and be sure the computer is communicating with the scanner. The PC will make a chime sound to verify communication has been established between the PC and scanner.



5. Go to SMI, Enter your Site ID



False Marcolle TD						
Enter Your Site ID:		Q				
Data Sent To Maximo Data N	OT Sant To Maximo	Data Fi	rom Scanner Unload			
Data Sent To Maximo Data N From Scanner Upload	OT Sent To Maximo ORGID	Data Fi	VA_AEMS_MERS_REFID	ASSETNUM	LOCATION	VA_LASTINVDATE
Data Sent To Maximo Data N From Scanner Upload	IOT Sent To Maximo ORGID	Data Fi SITEID	VA_AEMS_MERS_REFID	ASSETNUM	LOCATION	VA_LASTINVDAT

6. In the "Select Scanner Type" section, choose Other like shown in Figure 3.

		○ CS4070 ○ Other				
		om Scanner Upload	Data Fr	Sent To Maximo	Maximo Data NOT Ser	Data Sent To Maximo
LOCATION	ASSETNUM	VA_AEM5_MERS_REFID	SITEID	ORGID	r Upload	From Scanner Upload
	ASSETNUM	VA_AEM5_MERS_REFID	SITEID	ORGID	r Upload	From Scanner Upload

Figure 3

7. Choose Ready to Receive From Scanner

Data Sent To Maximo Data From Scanner Upload Data From Scanner Upload Location ValastinvDate ValastinVBY From Scanner Upload ORGID STEED ValastinVB ASSETNUM LOCATION ValastinvDate ValastinVBY	Data Sent To Maximo Data From Scanner Upload Data From Scanner Upload From Scanner Upload ORGID SITEID VA_AEHS_MERS_REFID ASSETNUM LOCATION VA_LASTINVDATE VA_LASTINVDATE VA_LASTINVDATE	Enter Your Site ID: NCHCS				Ready to Receiv From Scanner	re		
From Scanner Upload ORGID SITEID VA_AEHS_MERS_REFID ASSETNUM LOCATION VA_LASTINVDATE VA_LASTINVBY	From Scanner Upload ORGID SITEID VA_AEHS_MERS_REFID ASSETNUM LOCATION VA_LASTINVDATE VA_LASTINVDY	Data Sent To Maximo Data NOT Sent	To Maximo	Data From	1 Scanner Upload	ß			
		rom Scanner Upload	ORGID	SITEID	VA_AEMS_MERS_REFID	ASSETNUM	LOCATION	VA_LASTINVDATE	VA_LASTINVBY

8. Go to KWedge and click Transfer



9. Then choose Done Receiving in the SMI application

Enter Your Site ID: NCH	S				Do Rece	one eiving	
Data Sent To Maximo Data NOT From Scanner Upload	Sent To Maximo	Data Fr	VA_AEM5_MER5_REFID	ASSETNUM	LOCATION	VA_LASTINVDATE	VA_LASTINVBY
rom scanner opioau	OKGID	SITED	VA_AEPIS_PIEKS_KEPID	ASSEMUT	LOCATION	VA_CASIIMVDATE	VALUASTINVBT
			M				

Verify last line says "**END***. If it does not say END, then attempt another transfer. If END does appear, choose the SEND To Maximo to complete the process.

Enter Your Site ID: NCHCS	1	E.				Send	1
		1				To Max	imo
Data Sent To Maximo Data NOT Sent	To Maximo	Data Fr	rom Scanner Upload				
From Scanner Upload	ORGID	SITEID	VA_AEMS_MERS_REFID	ASSETNUM	LOCATION	VA_LASTINVDATE	VA_LASTINVBY
оміхам							
ID RRG							
SP103-R2-MTZ							
612 MX2557090,13:49:57,10/31/2022							
612 MX2557067,13:50:05,10/31/2022							
612 MX1911794,13:50:12,10/31/2022							
612 MX2631603,13:50:19,10/31/2022							
612 MX2557059,13:50:25,10/31/2022							
612 MX2557091,13:50:35,10/31/2022							

Recovering Scanned Data

Beginning with mobile app version 5.3.2, the scanners have begun retaining an archived copy of any Transfer file created to address concerns

Select Target	×
Current: Reflections Workspace - VISTA New target:	
KWedgel Microsoft Text Input Application Program Manager Reflections Workspace - VISTA	OK Cancel where

users believe scanned data has been lost. KWedge is not able to retrieve this file so the user must navigate to the file using Window's File Explorer and going to the folder path:

CT45 -> Internal shared storage -> MSS -> Export -> Archive

Archived files will be retained on the scanner until deleted. To delete one or multiple archived files, use a PC to go to the directory of the scanner and delete the files.



Special Messages

If you receive this message on your PC at any time, it means that the KWedge program does not have permission to delete the vista.txt.txt (vista.txt) file. Please have IT enable that capability as pointed out in Step 2 of the install instructions.

→Copy and Replace the fi	Replace le in the destination fold	er with the file you ar	e copying:
vista.txt.tx vista.txt (C:)	t (WINDOWS)		
Size: 656 by Date modifi	tes ied: 7/14/2020 11:55 AM		
→Don't copy No files will b	/ e changed. Leave this fil	e in the destination fo	lder:
vista.txt.tx vista.txt (C: Size: 656 by	t \MSS Software\KWedge ter	for Android)	
Data mandifi	indi 7/14/2020 11:55 AM		

Advanced Options

In some circumstances, KWedge can be configured to work differently than the standard setup. For example, the sound can be turned off so that it does not give you a "TA DA" when the transfer was successful. To get to this menu, select File | Options in KWedge. You will be presented with the following Options:

Options			?	×
Additional delay at start		seconds	Savea	and exit
Additional delay between characters	2	100ths of sec	Exit w/o	o saving
Additional delay between lines	5	100ths of sec		
Sounds Make sounds on critical err Make sounds on successful Make sounds on successful	ors ul devie	ce session		
 Use 12-hour notation for times Send status messages to log file, a 	also			
Show warnings when wedge is dis Suppress 'ENTER' on last record	abled	and a device co	onnects	

If you need help with configuring these options, please contact support@mss-software.com

Issues with Transferring Data

If you are having issues getting the data to transfer, the first thing to check is whether the scanner is communicating with a PC. The best way to check for a connection is to open your File Explorer and look under 'This PC'. The CK65 or CT45 will show up and if you click on it, you should be able to view its folders.

	~	Q	9	Search CT40
Internal shared storage 12.3 GB free of 17.8 GB	IPSM card 1.86 GB free of 1.90 GB			

Troubleshooting KWedge

- 1. Verify a Transfer File has been created. A common problem is forgetting to press option 2 from the Main Menu of the Scanner.
- 2. Verify the PC is communicating with the scanners. A common problem is an endpoint application such as McAfee is disallowing USB communication. Contact IT about whitelisting the scanner to allow USB communication.

Troubleshooting Checklist

- 1. Verify scanner is in its dock and the PC is communicating with it.
 - a. The PC should make a sound when the scanner is placed in its dock
 - b. You should be able to view the scanner's files using Windows File Explorer This PC -> CT45 (or CK65) -> Internal Shared storage

like graphic at right

2. Verify the inventory transfer file has been created as outlined in Step 4 of this document.



3. Verify the dock has power, the scanner's battery icon is showing it is charging

 Be sure KWedge for Android is being used on the PC, not KWedge 3.23 Be sure the correct target is selected. The target application is normally Reflections – VistA or Maximo.

Scan Count Pro aka SC Pro

SC Pro is a simple inventory collection program which creates an Excel file as the user scans and collects data. It may be configured in many ways such as with or without Location, with or without Quantity, and so forth. Each scanner can be configured uniquely from one another.

Situations where it may be helpful include scanning Model and Serial numbers of new goods being received in the warehouse

More Information and the User Manual may be found at the mss-software.com website here:



https://www.mss-software.com/shop/barcode_inventory_control_solutions/MSS-SCPRO-KIT1~Scan+My+Way+Inventory+System.htm



Home Page Layout Data File Formulas Revi View 🔏 Cut Calibri A A - 11 Copy ~ Paste ΙU H • 🔬 • ✓ Format Painter Clipboard N Font A2 > f_x L100 ÷ D B C F G A Bin 1 Item Qty 2 L100 987026 6 3 L100 980001 3 4 5 6 7 8

Sample Input Screen and Output File:

Call us at 800-428-8643 www.mss-software.com support@mss-software.com